



CERTIFICATE OF ANALYSIS



DETERGENT DISINFECTANT

Lot # 23920A

Date of Manufacture August 2020

Expiration Date August 2021

SIZE	16oz (475mL)	32oz (950mL)	64oz (1.9L)	1gal (3.8L)
DECON #	4101	4102	4103	4104

TEST	SPECIFICATION	UNIT	TEST RESULT
Active Content	0.21 - 0.23	% w/w	0.23
pH	10.7 - 12.7	N/A	12.1
Density	8.08 - 8.86 lb./gal	lb./gal	8.38
Appearance	Clear to Yellow Liquid Characteristic odor	N/A	Pass

Reviewed and Approved by:

Date: August 2020

A handwritten signature in blue ink, appearing to read 'Karen Wend'.

Karen Wend
Manager, Quality and Technical Services
Template revised 10/02/2012

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.