

Disposable SW Plate Riplate 384 - 200 µl



Product description

Riplate 384 - 200 µl


Part Number	43004-0200
Well volume	200 µl
Well bottom	V-Bottom
Plate color	clear
Sterile	no
Re-usable	no
Patient contact	no

Materials

Well PP – Polypropylene

Plate Dimensions

Dimension a	127.76 ± 0.25 mm
Dimension b	85.48 ± 0.25 mm
Dimension c	22.5 ± 0.2 mm
Dimension d	8.99 ± 0.13 mm
Dimension e	12.13 ± 0.13 mm
Dimension f	103.50 ± 0.25 mm
Dimension g	67.50 ± 0.25 mm
Dimension h	4.50 ± 0.08 mm
Dimension i	3.32 ± 0.13 mm
Dimension j	3.80 ± 0.13 mm
Dimension k	19.9 ± 0.2 mm
Dimension l	2.5 ± 0.1 mm



Packaging

Sealed bag

10 plates sealed

Riplate 384 - 200 µl

Shopping box (brown)

4 stacks (40 plates) in a box

REF 43004-0200 /	8°C / 2	
2025/10 2030/10	Menge: 40 Stück Quantity: 40 pieces	
LOT 1022592/1-1	1022592/1	64412872
	ritter MADE IN GERMANY	
Ritter GmbH - Kaufbeurer Str. 55 - 86830 Schwabmünchen - Germany		

Pallet

112 boxes (4480 plates)

Specifications are subject to change without prior notice.

In case of any questions, please contact our sales representative via phone, e-mail or our website.

<https://www.ritter-medical.de/en/contact>