

Product Specification

Test	Specification
Assay (by GC)	$\geq 95 \%$
Aromatic Compounds (as C ₆ H ₆)	$\leq 0.02 \%$
Color (APHA)	≤ 10
Density (g/ml) at 25°C	≤ 0.657
Residue after Evaporation	$\leq 0.001 \%$
Sulfur Compounds (as S)	$\leq 0.005 \%$
Thiophene	passes test
Water Soluble Titrable Acid (meq/g)	≤ 0.0003
Trace Impurities (in ppm)	
Aluminium (Al)	≤ 0.5
Barium (Ba)	≤ 0.1
Boron (B)	≤ 0.02
Cadmium (Cd)	≤ 0.05
Calcium (Ca)	≤ 0.5
Chromium (Cr)	≤ 0.02
Cobalt (Co)	≤ 0.02
Copper (Cu)	≤ 0.02
Iron (Fe)	≤ 0.1
Lead (Pb)	≤ 0.1
Magnesium (Mg)	≤ 0.1
Manganese (Mn)	≤ 0.02
Nickel (Ni)	≤ 0.02
Tin (Sn)	≤ 0.1
Zinc (Zn)	≤ 0.1

For Laboratory, Research or Manufacturing use.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000