



BOSTWICK CONSISTOMETER

RECOMMENDED FOR QUALITY ASSURANCE

Checks Laboratory or Production Samples
Against Consistency, Viscosity or
Flow Rate Standards

USES

■ Food Processing

- Tomato Sauces + Others
- Preserves
- Fillings
- Soups
- Baby Foods
- Salad Dressings
- Suggested by USDA

■ Paint Manufacturing

- Epoxy Removers
- Polyurethane Paint
- Adhesives
- Coatings
- Mil Spec R-81294B

■ Cosmetics Formulations

- Lotions • Shampoos
- Flow Base Creams

■ Chemical Production

- Selected Solvents
- Viscous Liquids
- Slurries
- Emulsions



OPERATIONAL FEATURES & SPECIFICATIONS

■ Stainless Steel Construction

- Resists Corrosion
- Holds Up Under Plant Conditions
- Appropriate for Laboratory Use

■ Engraved Graduations In 0.5 cm Divisions

- Wear & Smear Resistant
- Assures Accurate Results

■ Small Size 14" L x 3½" W x 5½" H

- Uses Little Bench Space
- Requires 75 ML Sample

■ Spring Loaded Gate

- Instantaneous Gate Release
- Prevents Premature Flowing of Sample

■ 2 Leveling Screws & Spirit Level

- Precise Positioning of Flow Angle
- Assures Reproducibility



CSC SCIENTIFIC COMPANY, INC.

MANUFACTURERS AND DISTRIBUTORS OF LABORATORY EQUIPMENT