

Reagecon

Delivering the Correct Result



Reagecon Standard Colour Solution Y (Yellow) according to European Pharmacopoeia (EP) Chapter 2.2.2

Item No. EP706

Reagecon's Standard Colour Solution Y (Yellow) is formulated as a reference solution for the examination of the degree of colouration of liquids in accordance with the European Pharmacopoeia Chapter 2.2.2. This product is produced by combining Reagecon's Primary Colour Solutions.

Additional Details

Pack Size	100ml
Weight	0.15
Weight Unit	Kg
Colour	Yellow
Tariff Code	3822 190090
Country of Origin	Ireland

Safety Information

UN Number	1760
Class	8
GHS	

Price € 110,¹¹ price per 100ml
excluding tax

Additional Details

Pack Size	100ml
Weight	0.15
Weight Unit	Kg
Colour	Yellow
Tariff Code	3822 190090
Country of Origin	Ireland

GHS

DESCRIPTION

This Standard Colour Solution is a ready to use coloration standard presented in high quality, tamper evident bottles for use as a calibration and/or quality control standard. It is prepared by combining Reagecon Primary Colour Solutions.

Manufactured gravimetrically according to European Pharmacopoeia Chapter 2.2.2 guidelines with a Balance Certificate of Calibration issued in accordance with ISO 17025 (INAB ref: 265C)

Concentration verified using a high performance spectrophotometer.

PRODUCT SPECIFICATIONS

European Pharmacopoeia (EP) Standard, European Pharmacopoeia (EP) Material

Use	Colouration Standard	Packaging	Amber Glass Winchester Bottle
Tested	Yes	Physical Form	Liquid