



Certificate of Conformance, Analysis, and Sterility

This certification is provided as full assurance that the following product code and lot number were manufactured in accordance via a Quality System certified to ISO 9001 with prescribed procedures and specifications.

Product Code:	PS-911EB/ETOH
Lot Number:	1208341-0-0-1
Manufacturing Date:	02/02/2026
Ethyl Alcohol(denatured):	70 ± 3% by volume with the remainder DI H2O.
Expiration Date:	02/2028

Analysis

Quantity of Particles Po:	1.127 millions/m ²
---------------------------	-------------------------------

Sterility

Gamma Processing Run ID:	3915882
Irradiation Dose Specified:	18.2 kGy to 40.0 kGy
Irradiation Dose Delivered:	19.0 kGy to 33.6 kGy
Validated Sterile Product:	Validation Sterile packets available upon request.
Sterility Assurance Level:	10 ⁻⁶ (SAL determined by test method AAMI/ISO 11137)

QC Supervisor
coc@contecinc.com

Monday, February 16, 2026

NORTH AMERICA

Contec, Inc.
525 Locust Grove
Spartanburg, SC 29303 USA
tel: +1 864 503 8333
fax: +1 864 503 8444

EUROPE

Contec Europe
ZI Du Prat, Avenue Paul Dupleix
56000 VANNES France
tel: +33 (0)2 97 43 76 98
fax: +33 (0) 97 43 76 86

ASIA

Contec Cleanroom Technology (Suzhou) Co., Ltd.
17 Longyun Road
Suzhou Industrial Park, Suzhou 215024 China
tel: +86 512 6274 4050
fax: +86 512 6274 4051

web: www.contecinc.com
e-mail: info@contecinc.com

Page 1 of 1

DID: 2061B



Certificate of Processing

STERIGENICS 10811 Withers Cove Park Dr. Charlotte NC 28278
 TEL 704 588-6877 FAX 704 588-8355 www.sterigenics.com

R55480102

02/11/2026 13:19:55 GMT
 Page 1 of 1

Customer Name:	Contec- PO Box 530	Processing Facility:	Charlotte Gamma	Work Order #	3915882
P.O.#	156401			Sales Order #	3997537
	18.2-40.0 kGy	CNT4, Gamma Treatment	Irradiation Date/Time:	02/10/2026 16:57:00 GMT	

SO Line #	Qty	UOM	Customer Item Number	Customer Item Description	Customer Lot Number	Customer Load Number
101.000	414	CA	PS-911EB/ETOH	Clean Room Product	1208341-0-0-1	156401
	414	CA	Total			

Quality Test Summary

Op#	Quality Test Description	Minimum Spec	Maximum Spec	Result	Pass/Fail	-----Signed By----- User	Date /Time
450.00	Minimum Dose	18.2 kGy	40.0 kGy	19.0 KGY	Pass	JALCIVAR	02/11/2026 03:41:35 GMT
		Reason Code Test				JOHANN ALCIVAR	
450.00	Maximum Dose	18.2 kGy	40.0 kGy	33.6 KGY	Pass	JALCIVAR	02/11/2026 03:41:46 GMT
		Reason Code Test				JOHANN ALCIVAR	

Sterigenics certifies that the materials listed above (as described by the Manufacturer) received the indicated doses within the precision and accuracy of the dosimetry system employed.

Electronically Signed By: TAMMY MCDANIEL
 Reason: Work Order Completion

Date: 02/11/2026 13:17:29 GMT