

TECHNICAL DATA SHEET

DESCRIPTION - Luer Filling Cannula / Needles



REF – Issue 1

REF; Various

MASTER SAMPLE YES

PRODUCTION SAMPLE

COMPONENT PART	QTY	MATERIAL
Mount	1	¼" Sq Brass CZ121 No 108 or 109 mount UK Supply
Cannula / Needle tube	1	Stainless Steel tubing type 304 Japan Supply

FINAL FINISH

POLISHED X

NICKEL X

CHROME

SILVER

INSPECTION CRITICAL POINTS

Rounded smooth face on cannula tube
Burr free point on needle tip
Debris and burr free bore on both styles

CRITICAL SAFETY FEATURES

VARIANTS

DETAILS

No 108 mount

Small mount used for gauges up to 20g

No 109 mount

Large mount used for gauge 19 and upwards

PACKAGING

Covered by white plastic straw (cut to length)

Packed in multiples of 10 in 3" lay flat tubing 500g