



Product Specification Sheet

VWR STORAGE RACKS FOR SLIDE BOXES

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Cat. No.	:	631-0742
Material	:	Acrylic
Power	:	N/A
Colour	:	Clear
Number Held	:	15 Slide Boxes
Array	:	3 Columns x 5 Rows
Store	:	Can fit either 1 x 100 place slide box, or 2 x 50 place slide boxes per slot in this Acrylic Storage Rack, so the maximum slide capacity for the complete rack is 1500 slides using either 50 or 100 place, or a combination of both types of slide boxes.
Dimensions	:	711 x 171 x 28 mm
Slot Depth	:	166 mm
Standards	:	N/A

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.

Date of Issue: **May, 2011**