

[Bowls & Dishes \(index.php?gruppen&kategorie=0010&Seite=1\)](index.php?gruppen&kategorie=0010&Seite=1)

[Miscellaneous Items \(index.php?gruppen&kategorie=0011&Seite=1\)](index.php?gruppen&kategorie=0011&Seite=1)

[Labjacks \(index.php?gruppen&kategorie=0012&Seite=1\)](index.php?gruppen&kategorie=0012&Seite=1)

[Water Jet Pumps \(index.php?gruppen&kategorie=0013&Seite=1\)](index.php?gruppen&kategorie=0013&Seite=1)

## Specifications

| Item no. | temperature<br>°C | consumption<br>l/h | height<br>mm | weight<br>g | tube<br>Ø mm | head<br>Ø mm | load<br>kW | DIN | GS<br> | Natural gas | Bottle gas | safety-base<br>+ 38 g |
|----------|-------------------|--------------------|--------------|-------------|--------------|--------------|------------|-----|-------------------------------------------------------------------------------------------|-------------|------------|-----------------------|
| 1000     | 1500              | 110                | 155          | 210         | 14           | 17           | 1,25       | +   | +                                                                                         | +           |            | •                     |
| 1001     | 1540              | 55                 | 155          | 210         | 14           | 17           | 1,35       | +   | +                                                                                         |             | +          | •                     |
| 1010     | 1500              | 110                | 155          | 220         | 14           | 17           | 1,00       |     |                                                                                           | +           |            | •                     |
| 1010/1   | 1500              | 115                | 155          | 220         | 14           | 17           | 1,00       | +   | +                                                                                         | +           |            | •                     |
| 1011     | 1550              | 55                 | 155          | 220         | 14           | 17           | 1,70       |     |                                                                                           |             | +          | •                     |
| 1011/1   | 1550              | 55                 | 155          | 220         | 14           | 17           | 1,70       | +   | +                                                                                         |             | +          | •                     |
| 1020     | 1500              | 110                | 155          | 220         | 14           | 17           | 1,10       | +   | +                                                                                         | +           |            | •                     |
| 1021     | 1520              | 55                 | 155          | 220         | 14           | 17           | 1,25       | +   | +                                                                                         |             | +          | •                     |
| 1030     | 1500              | 110                | 155          | 225         | 14           | 17           | 1,10       |     |                                                                                           | +           |            | •                     |
| 1030/1   | 1500              | 110                | 155          | 225         | 14           | 17           | 1,10       | +   | +                                                                                         | +           |            | •                     |
| 1031     | 1520              | 50                 | 155          | 225         | 14           | 17           | 1,25       |     |                                                                                           |             | +          | •                     |
| 1031/1   | 1520              | 50                 | 155          | 225         | 14           | 17           | 1,25       | +   | +                                                                                         |             | +          | •                     |
| 1040     | 1500              | 100                | 155          | 245         | 14           | 17           | 1,10       |     |                                                                                           | +           |            | •                     |
| 1040/1   | 1500              | 100                | 155          | 245         | 14           | 17           | 1,10       | +   | +                                                                                         | +           |            | •                     |
| 1041     | 1520              | 55                 | 155          | 245         | 14           | 17           | 1,25       |     |                                                                                           |             | +          | •                     |
| 1041/1   | 1520              | 55                 | 155          | 245         | 14           | 17           | 1,25       | +   | +                                                                                         |             | +          | •                     |
| 1045/1   | 1480              | 90                 | 155          | 235         | 14           | 17           | 1,10       | +   | +                                                                                         | +           |            | •                     |
| 1046/1   | 1500              | 55                 | 155          | 235         | 14           | 17           | 1,25       | +   | +                                                                                         |             | +          | •                     |
| 1050     | 1460              | 110                | 170          | 270         | 13           | 17           | 1,20       |     |                                                                                           | +           |            | •                     |
| 1050     | 1510              | 35                 | 170          | 270         | 13           | 17           | 1,20       |     |                                                                                           |             | +          | •                     |
| 1060     | 1600              | 120                | 125          | 120         | 14           | 17           | 1,30       |     |                                                                                           | +           |            | •                     |
| 1061     | 1630              | 42                 | 125          | 120         | 14           | 17           | 1,00       |     |                                                                                           |             | +          | •                     |
| 1110     | 1310              | 35                 | 100          | 115         | 7            | 11           | 0,40       |     |                                                                                           | +           |            | •                     |
| 1111     | 1320              | 20                 | 100          | 115         | 7            | 11           | 0,30       |     |                                                                                           |             | +          | •                     |
| 1200     | 1550              | 115                | 155          | 245         | 14           | 17           | 1,50       |     |                                                                                           | +           |            | •                     |
| 1201     | 1560              | 65                 | 155          | 245         | 14           | 17           | 1,30       |     |                                                                                           |             | +          | •                     |
| 1210     | 1520              | 125                | 155          | 260         | 14           | 17           | 1,30       |     |                                                                                           | +           |            | •                     |
| 1210/1   | 1520              | 125                | 155          | 260         | 14           | 17           | 1,30       | +   | +                                                                                         | +           |            | •                     |
| 1211     | 1560              | 65                 | 155          | 260         | 14           | 17           | 1,50       |     |                                                                                           |             | +          | •                     |
| 1211/1   | 1560              | 65                 | 155          | 260         | 14           | 17           | 1,50       | +   | +                                                                                         |             | +          | •                     |
| 1220     | 1530              | 110                | 155          | 260         | 14           | 17           | 1,25       | +   | +                                                                                         | +           |            | •                     |
| 1221     | 1530              | 65                 | 155          | 260         | 14           | 17           | 1,50       | +   | +                                                                                         |             | +          | •                     |