

VWR® BACTI CELL SPREADERS

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	60828-682
Manufacture Lot Number:	021320
Sterilization Lot Number:	17645
Date of Sterilization:	March 16, 2020
Dose Range Received:	25-25 kGy
Description:	VWR Spreader 30mm Pk25.
Country of Origin:	USA

Quality System Compliance

Products are manufactured under GMP guideline standards. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

Product Specifications

Class I Device

Material: Virgin Polypropylene

Sterilization Process : E-Beam

BSE/TSE: The product does not contain or was manufactured with any animal products or materials of animal origin.

RNase/DNase/PYROGEN: The product is RNase/DNase/PYROGEN free.

Cytotoxicity: The product is non cytotoxic.

Latex Free: The product is latex free.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 16, 2020