

## CERTIFICATE OF CONFORMITY & CERTIFICATE OF ANALYSIS

|                    |                                                                            |
|--------------------|----------------------------------------------------------------------------|
| Brand              | MICROFLEX®                                                                 |
| Size(s)            | XS (5.5 - 6),S (6.5 - 7),M (7.5 - 8)                                       |
| Material           | Nitrile                                                                    |
| Style              | XCEED® XC-310                                                              |
| Short Description  | Thin, Strong Nitrile Exam Glove with ERGOFORM™ Ergonomic Design Technology |
| Lot Number         | 2501PM                                                                     |
| Batch Number       | 2501037IPM                                                                 |
| Manufacturing Date | 2025-01                                                                    |
| Expiration Date    | 2030-01                                                                    |

### CERTIFICATE OF CONFORMITY

Ansell certifies that the above product is manufactured in accordance with Ansell product specifications and is in conformity with the following standards or normative documents:

| Standard/Normative Document | Result   |
|-----------------------------|----------|
| FDA 21 CFR 177-2600         | Conforms |
| Reach Regulation            | Conforms |

### CERTIFICATE OF ANALYSIS

| Property                                               | Test Method           | Specification | Result   |
|--------------------------------------------------------|-----------------------|---------------|----------|
| <b>Freedom From Holes (General Inspection Level I)</b> |                       |               |          |
| Freedom from Holes - Final Release                     | EN455-1/ASTM D5151    | G1, AQL 0.65  | Conforms |
| <b>Visual Inspection</b>                               |                       |               |          |
| Visual Inspection-major                                | Internal test method  | G-I, AQL 2.5  | Conforms |
| Visual Inspection-minor                                | Internal test method  | G1, AQL 4.0   | Conforms |
| <b>Physical Properties</b>                             |                       |               |          |
| Material Standard                                      | ASTM D6319            | S2, AQL 4.0   | Conforms |
| Dimension Test                                         | ASTM D6319/ASTM D3767 | S2, AQL 4.0   | Conforms |
| Finger Thickness (mm)                                  | ASTM D6319/ASTM D3767 | 0.09 - 0.13   | Conforms |
| Palm Thickness (mm)                                    | ASTM D6319/ASTM D3767 | 0.05- 0.09    | Conforms |
| Cuff Thickness (mm)                                    | ASTM D6319/ASTM D3767 | 0.04 -0.08    | Conforms |
| Glove Length (mm)                                      | ASTM D6319/ASTM D3767 | Min. 240      | Conforms |
| Palm Width (mm)                                        | ASTM D6319/ASTM D3767 | 71-120        | Conforms |
| Physical Properties                                    | ASTM D412             | S2, AQL 4.0   | Conforms |
| Tensile Strength (MPa)                                 | ASTM D412             | Min. 14       | Conforms |
| Ultimate Elongation (%)                                | ASTM D412             | Min. 400      | Conforms |
| Modulus at 300% elongation (MPa)                       | ASTM D412             | Max. 10       | Conforms |

| Property                           | Test Method | Specification | Result   |
|------------------------------------|-------------|---------------|----------|
| <b>Additional Properties</b>       |             |               |          |
| Residual Powder Content (mg/glove) | ASTM D6124  | Max. 2.0      | Conforms |
| Protein Content (µg/g)             | ASTM D5712  | N/A           | Conforms |

Electronically Signed by:

**Sandra Wessel-Jess**

QA Manager - Americas

Reno, NV, United States

Date: 23 October 2025

For additional information visit us at [www.ansell.com](http://www.ansell.com), or contact us at:

**Europe, Middle East & Africa Region**

Ansell Healthcare Europe N.V.

Riverside Business Park

Blvd International, 55

B-1070 Brussels,

Belgium

T: +32 (0) 2 528 74 00

F: +32 (0) 2 528 74 01

**Australia & Asia Pacific Region**

678 Victoria St,

Richmond VIC 3121

Australia

Tel: +61 1800 337 041

Fax: +611 800 803578

Email: [anzcustomerservice@ansell.com](mailto:anzcustomerservice@ansell.com)

Prima 6, Prima Avenue, Block 3512

Jalan Teknokrat 6

63000 Cyberjaya, Malaysia

Tel: +603-8310 6688

Fax: +603-8310 6699

Email: [cyberjaya@ansell.com](mailto:cyberjaya@ansell.com)

**North America Region**

Ansell Healthcare products LLC

111 S Wood Ave,

Suite 900 Iselin,

NJ 08830,

United States

T: +1 800 800 0444

F: +1 800 800 0445

CA Tel: +1 800 363 8340

Email: [customerserviceus@ansell.com](mailto:customerserviceus@ansell.com)

Ansell Healthcare Products LLC

2301 Robb Drive,

Reno,

NV 89523,

United States

US Tel: +1 800 800 0444

US Fax: +1 800 800 0445

Email: [customerserviceus@ansell.com](mailto:customerserviceus@ansell.com)

**Latin America & Caribbean Region**

Ansell Comercial México S.A. de C.V.

Blvd. Bernardo Quintana 7001-C,

Q7001 Torre 2, Suites 1304-1306 Centro Sur.

Querétaro, Qro. México C.P. 76090

Tel: +52 442 296 2050

Email: [cslac@ansell.com](mailto:cslac@ansell.com)

**Ansell Brazil Ltda.**

Rua das Figueiras, 474

4º andar - Bairro Jardim

Santo André - SP

CEP: 09080-300 – Brasil

Tel: +55 11 3356 3100

Email: [sac@ansell.com](mailto:sac@ansell.com)