

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Test tubes with or without rim, thick walled
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2023-03-16

With Rim

Catalogue Number	Length [mm]	Ø ext [mm]	Wall Thickness [mm]
212-0766	75	12	0.9 – 1.0
212-0767	100	12	0.9 – 1.0
212-0768	130	14	0.9 – 1.0
212-0769	130	16	1.1 – 1.2
212-0772	150	20	1.1 – 1.2
212-0774	150	25	1.1 – 1.2
212-0770	160	16	1.1 – 1.2
212-0771	180	18	1.1 – 1.2
212-0773	180	20	1.1 – 1.2
212-0775	200	25	1.1 – 1.2
212-0776	200	30	1.1 – 1.2

Without Rim

Catalogue Number	Length [mm]	Ø ext [mm]	Wall Thickness [mm]
212-0777	75	10	0.8 – 1.0
212-0778	75	12	0.8 – 1.0
212-0779	100	10	0.8 – 1.0
212-0780	100	12	0.8 – 1.0
212-0784	130	16	1.1 – 1.2
212-0783	130	14	0.9 – 1.0

Catalogue Number	Length [mm]	Ø ext [mm]	Wall Thickness [mm]
212-0785	150	20	1.1 – 1.2
212-0787	150	25	1.1 – 1.2
212-0781	160	16	0.8 – 1.0
212-0782	180	18	1.1 – 1.2
212-0786	180	20	1.1 – 1.2
212-0788	200	25	1.1 – 1.2
212-0789	200	30	1.1 – 1.2

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard and calibrated in **ISO 17025:2017** certified Laboratory. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Technical standards for the products: NA

Calibration procedure: NA

Material	The products are made from glass tubing that conforms to ASTM E-438 Type-1, Class A Borosilicate Glass 3.3 Expansion.
Sterilization	Hot air sterilization, in the oven up to 140 °C Steam sterilization, in an autoclave 121 °C/ 20 min 134 °C/ 10 min.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting

Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	03-2023
2	1.1	Additional Product code included	06-2023