

VWR® 600mL Treated Tissue Culture Flasks

Description: 600ML, Treated, w/plug seal cap, Sterilized
600ML, Treated, w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired

Materials



Flask: GPPS (General polystyrene) Color: Clear

Cap: HDPE(High Density Polyethylene) Color: Purple 259

Membrane: 0.22µm Hydrophobic PP (Only use for Vented Cap)


Features

- Temperature: -20°C~+50°C
- Sterilation: sterilized by gamma irradiation SAL10⁶(ISO11137)
- Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Working volume: 400 ml
- Growth area: 182.0 cm²
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production



VWR® 600ml Treated Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2316	10062-866	600ML Tissue Culture Flask, Treated, w/plug seal cap, Sterilized	Yes	5	20	40
734-2315	10062-864	600ML Tissue Culture Flask, Treated, w/vent cap, Sterilized	Yes	5	20	40



VWR® 850mL Treated Tissue Culture Flasks

Description: 850ML, Treated, w/plug seal cap, Sterilized
850ML, Treated, w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is desired



Materials

- Flask: GPPS (General polystyrene) Color: Clear
- Cap: HDPE(High Density Polyethylene) Color: Purple 259
- Membrane: 0.22µm Hydrophobic PP (Only use for Vented Cap)

Features

- Temperature: -20°C~+50°C
- Sterilation: sterilized by gamma irradiation SAL10⁶(ISO11137)
- Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Working volume: 800 ml
- Growth area: 300.0 cm2
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production

Dimensions (Unit: mm)



VWR® 850ml Treated Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2601	10062-886	850ML Tissue Culture Flask, Treated, w/plug seal cap, Sterilized	Yes	3	9	18
734-2600	10062-884	850ML Tissue Culture Flask, Treated, w/vent cap, Sterilized	Yes	3	9	18