

Certificate of Conformity

VWR Catalogue Number	VWRA112-5399
Description	VWR®, Latex Gloves, Powder-free, Non-Sterile, M, 240 mm length, 65 mm thickness, 122 mm diameter
Lot Number	K01251455
Date of Manufacture (yyyy-mm-dd)	2025-01
Expiry Date	2027-12
Country of Origin	Manufactured in Malaysia
Date of Issue (yyyy-mm-dd)	2026-06-16

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Ambidextrous / hand-specific	Ambidextrous
Disposable / reusable	Disposable
Powdered / powder-free	Powder-free
Latex	Latex-free
Cuff style	Beaded
Cuff thickness	0.10 mm*
Finish	Textured

Material Natural rubber latex.

Sterilization Non-sterile.

ATP Assay Not Applicable.

VWR Singapore Pte Ltd, 18 Gul Drive, Singapore 629468 - <http://www.vwr.com>

DNase & RNase Free	Not Applicable.
BSE/TSE	Products does not contain tissues, cells and their derivatives of animal origin (Regulation (EU) No. 722/2012).
Cytotoxicity	Product showed cytotoxic reactivity at the 0.4 cm ² /mL extract concentrations and no cytotoxic reactivity at the 0.2 cm ² /mL extract concentrations under the condition of this test. Cytotoxicity results acceptable intended use (short-term use on intact skin).
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Product contains natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP has not been intentionally incorporated in the manufacturing of the products.
RoHS	Not Applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.
Storage Conditions	Store in original packaging in a dry and dark place at 10° C to 30° C. Refer to guidelines of storage of rubber products as described in ISO 2230:2002. Ensure that storage area is kept cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper ions discolour the glove. Protect gloves against ultraviolet light sources, such as sunlight and oxidizing agents. Storage above 30° C will lead to accelerated aging and should be avoided.
Shelf life	3 Years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Annexure

Parameter	Specification	Result	Sampling size
1.Dimension			
Palm width	95 ± 5	Min. 94 Max: 96	Pass
Length	Median : 230 mm	Median : 244 mm	Pass
Cuff Thickness	Min. 0.06 mm	Min. 0.09 mm	Pass
Palm Thickness	Min. 0.10 mm	Min. 0.12 mm	Pass
Finger Thickness	Min. 0.11 mm	Min. 0.14 mm	Pass
2. Functionality and Appearance test			
Barrier	G1 AQL 1.5	Pass	
Visual	G1 AQL 2.5	Pass	

Parameter	Specification	Result	Sampling size
3. Lab testing			
Tensile Strength (Un-aged) (Mpa)	Min : 18 Mpa	20.9	Pass
Tensile Strength (Aged) (Mpa)	Min : 14 Mpa	15.7	Pass
Elongation at Break (Un-aged) (%)	Min : 650%	720	Pass
Elongation at Break (Aged) (%)	Min : 500%	652	Pass
Powder (mg)	≤ 2 mg	0.74	Pass