

Material Number:

85415.601

Material description:

pH test strips 7.0-14.0 dosatest

Date specification set:15th Feb 2018

(All batches of this product after this date are tested to this specification. In the unlikely event that the specification is changed a new certificate will be issued)

Description

The colour fixed indicator strips allow for the determination of pH values between 7.0 and 14.0 in aqueous solutions. The indicator strip contains 3 test pads.

Presentation

Content: 100 indicator strips

Package dimensions: 66 x 90 x 10 mm

Indicator strip dimensions: 5.5 x 85 mm

Test pad dimensions: 5.0 x 5.5 mm.

Gradations

pH Values: 7.0 - 7.5 - 8.0 - 8.5 - 9.0 - 9.5 - 10.0 - 10.5 - 11.0 - 11.5 - 12.0 - 12.5 - 13.0 - 13.5 – 14.0.

Shelf Life

4 years.

Storage Temperature

4-30 ° C.

Storage Conditions

Avoid exposing the test strips to sunlight and moisture. Keep container cool and dry.

Signed on behalf of VWR:

Ben Petri

Quality Manager