



Specification

1.07043.1000 Universal peptone M 66 for microbiology

Specification		
Peptone content (ex N, calc. on dried substance)	≥ 80	%
Amino nitrogen (as N)	≥ 3.5	%
pH-value (2 %; water)	6 - 7	
Sulfated ash (800 °C)	≤ 15	%
Heavy metals (as Pb)	≤ 0.001	%
Loss on Drying (105 °C)	≤ 5	%
Suitability for microbiology	passes test	

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.