

Certificate of Conformity

VWR Catalogue Number	See List
Description	Temperature indicating strips
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2021-06-11

With ± 1.0 °C Accuracy

Catalogue Number	Range [°C]
620-9104	116 to 154
620-9105	160 to 199
620-9106	204 to 260
620-9102	37 to 65
620-9103	71 to 110

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Strips are printed laminates composed of self-adhesive films (with modified acrylic adhesives) polyester film or polyimide film (according to type), paper and various organic crystalline waxes.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Product have been shown to be non-toxic.
Non-Pyrogenic Statement	Product is free from Endotoxin.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store products containing these temperatures sealed in a plastic bag in a refrigerator ($\leq 10^{\circ}\text{C}/50^{\circ}\text{F}$).
Shelf Life	12 months.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.