

5 Technical data

Dimensions (L x W x H)	420x306x34 mm
Heating plate	CERAN-ceramic; heated area with grid, 200x200 mm
Power voltage	115/230 V, +/- 10%; 50/60 Hz
Max. Power consumption	400 W
Temperature control range	25 – 200°C
Fuses	315 mA, 250 V slow blow.

6 Temperature curves




Fig. 2 Temperature curves for heating and cooling

CAMAG (Switzerland) · Sonnenmattstrasse 11 · CH-4132 Muttenz 1

Telephone +41 61 467 34 34 · Fax +41 61 461 07 02 · E-Mail: info@camag.com

CAMAG (Germany) · Bismarckstraße 27-29 · D-12169 Berlin

Telephone +49 30 516 55 50 · Fax +49 30 795 70 73 · E-Mail: infoberlin@camag.com

CAMAG Scientific Inc. (USA) · 515 Cornelius Harnett Drive · Wilmington, NC 28401

Telephone 800 334 3909 · Fax 910 343 1834 · E-Mail: tlc@camag.com

www.camag.com

EC – Declaration of Conformity

We, CAMAG Chemie-Erzeugnisse und Adsorptionstechnik AG
Sonnenmattstrasse 11
4132 Muttenz
Switzerland

declare under our sole responsibility that the product

CAMAG® TLC Plate Heater 3

Product name

022.3306/ 022.3307

Article number(s)

to which this declaration relates is in conformity with the following provisions of directive(s):

- 2006/95/EC
- 2004/108/EC

Following standard(s) or other normative document(s):

- EN61010-1: 2010
- EN61326-1: 2013

Muttenz, 26 March 2015



Walter Rahm, Head of Quality Management

SWISS
MADE 

CE

UKCA – Declaration of Conformity

We, CAMAG Chemie-Erzeugnisse und Adsorptionstechnik AG
Sonnenmattstrasse 11
4132 Muttenz
Switzerland

declare under our sole responsibility that the product

CAMAG® TLC Plate Heater

Product name

022.3306 / 022.3307

Article number(s)

to which this declaration relates is in conformity with the following provisions of directive(s):

2014/35/EU LVD
2014/30/EU EMC
2011/65/EU RoHS

Following standard(s) or other normative document(s):

IEC 61010-1:2010/A1:2019/AC:2019-04
EN IEC 61010-2-010:2020
IEC 61326-1:2013

Muttenz, 12 January 2023

V. Maire-Widmer

Valeria Maire-Widmer, Head Quality Management

