

Residual Metal Catalyst/Reagents - Declaration

20904 2-PROPANOL PH.EUR./USP

With reference to guideline EMEA/CHMP/SWP/4446/2000, we hereby confirm that:

Pt, Pd, Ir, Rh, Ru, Os are below 0,5 ppm (Group) and Mo, V, Ni, Cr are below 1,0 ppm (Each)

Cu, Mn are below 1,5 ppm (Each)

Zn, Fe are below 2,0 ppm (Each)



Evelyn Verstraete
QA Private Label Chemicals



Esther Venmans
Manager QA
Chemical Manufacturing Europe

24/07/2017