



BRAND GMBH + CO KG

114832 BRAND Cryogenic Tubes
with external thread

Material:	Polypropylen, USP class VI and FDA compliant
Volume:	2.0 ml
Graduation up to:	1.8 ml
Sterility:	sterile, free from detectable pyrogens/endotoxins, DNA, DNase/ RNase and cytotoxic substances
Type:	self-standing
Tube diameter:	12,5 mm
Temperature stability:	+121°C to -196°C
Storage conditions:	originally packed, dry, protected from direct sun light
Shelf life:	60 month from day of manufacturing
CE marked (MD, IVD):	no
Batch controlled:	yes

Dimensions (in mm):

