

Certificate of Conformity

VWR Catalogue Number	See List
Description	Biohazard Bags, 50 µm Thickness
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-09-27

Red, with printed sign and temperature indicator, DOT approved

Catalogue Number	Capacity [L]	Dimension [mm]
129-1555	-	305 × 200
129-0033	12	610 × 305
129-1556	27	610 × 483
129-1557	-	762 × 610
129-0035	75	914 × 610
129-0034	80	890 × 635
129-1558	95	914 × 762
129-0036	260	1219 × 965

Plain red without printed sign

Catalogue Number	Capacity [L]	Dimension [mm]
129-0037	12	610 × 305
129-1559	49	762 × 610
129-0038	75	914 × 610

Biohazard red

Catalogue Number	Capacity [L]	Dimension [mm]
129-1550	-	305 × 200
129-1551	11	480 × 355
129-1552	27	610 × 483
129-1553	80	889 × 635

Catalogue Number	Capacity [L]	Dimension [mm]
129-1554	200	1143 × 914

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Established physical properties are tested in accordance with internationally recognized ASTM Test methods for Dart Impact Strength (ASTM D1709) and tear resistance (ASTM D1922).

Product Specifications

Material	Polypropylene Copolymer. Supplied with 12.7 cm twist ties. Thickness of bags is 50 µm.
Sterilization	Not Applicable. Products can be autoclaved at operating temperatures of 275 °F (135 °C).
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Resins and films used to manufacture this product may contain byproducts of bovine sourced raw material in their manufacture. These animal byproduct additives are processed under highly refined and controlled industrial chemical processes and at very high temperatures prior to purification and subsequent use as a polymer additive. The animal derived byproducts are not added during the process of manufacturing, but originate either in the manufacture of the resins or in the additives used in the film formulations.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Please see RoHS information below.

RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Room temperature.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.