

Technical data sheet

Artikel/article: Burette Class B, black graduation

Artikel-Nr./code no.: 310 01 60 60

Material/material: Borosilicate glass 3.3, DURAN

Class: B

Inhalt/ content: 10 ml

Tolerance: $\pm 0,03$ ml

Verpackungseinheit/
Packing unit: 2 pcs.

Standard: ISO 385

Hirschmann Laborgeräte GmbH & Co.KG
Eberstadt, 23.01.2024



Quality Manager, 31.01.2024

