



Funnels

Analytical & Disposable Funnel



Powder Funnel



Funnel Product Details

VWR Cat #	EU Cat#	Type of Funnel	Size	Material	Color	Funnel Outer Dia	Stem Dia	Stem Height	Total Height
414004-269	221-0239	Pow der	65 mm	PP	Natural	70	15.6	20.7	66.7
414004-270	221-0240	Pow der	80 mm	PP	Natural	84.8	15.8	30	85.2
414004-271	221-0241	Pow der	100 mm	PP	Natural	106.5	20.5	33	106.2
414004-272	221-0242	Pow der	150 mm	PP	Natural	149.4	26.7	33	136.5
414004-273	221-0243	Analytical	35 mm	PP	Natural	41.1	6.3	49.3	75.8
414004-274	221-0244	Analytical	45 mm	PP	Natural	50.3	7.1	52	84.6
414004-275	221-0245	Analytical	55 mm	PP	Natural	61.5	7.1	59.5	103.4
414004-276	221-0246	Analytical	65 mm	PP	Natural	71.1	6.4	61.4	111.7
414004-277	221-0247	Analytical	75 mm	PP	Natural	83.2	6.8	76.3	139.4
414004-278	221-0248	Analytical	90 mm	PP	Natural	98.5	8.2	89	164
414004-279	221-0249	Analytical	100 mm	PP	Natural	108.7	9.1	95.4	180
414004-280	221-0250	Analytical	160 mm	PP	Natural	168	14	154	280
414004-288	221-0251	Disposable	65 mm	PP	Natural	68	7.4	22.7	74.4
414004-289	221-0252	Disposable	55 mm	PP	Natural	60	7.4	27.1	71.3
414004-290	221-0253	Disposable	55 mm	PS	Natural	60	7.4	27.1	71.3

Material Recycle Code



Description: Funnels

Drawing:
Not to Scale



Release: 9/5/2011

Revision: 0