



Certificate of Analysis

Transferrin, Holo, Human Plasma, Low Endotoxin

Batch Number: 3534067
Material Number: 616424-100MG
Molecular Weight: 80000
CAS Number: 11096-37-0

Quality Release Date: 29 OCT 2020
Recommended Retest Date: 16 MAY 2022

Analytical Data

Test	Tolerance	Result
Formulation:		Lyophilized in 20 mM sodium phosphate, pH 7.4, 150 mM NaCl
Solubility:		H ₂ O. Please refer to the vial label for the pre-lyophilization concentration.
Form:		Powder
SDS PAGE purity:	≥95.0 %	≥95.0 %
"HBsAg; HIV,HCV antibodies":		Negative
Concentration before lyophilization:		120 mg/mL
Endotoxins:	≤0.1 EnU/mg	<0.1 EnU/mg
Storage and Handling:	-20°C	

This lot conforms to specifications established by EMD Millipore Corporation for this product.

Quality Control/ Assurance Signature

29 OCT 2020
Date

Prices and availability are subject to change. ©Copyright 2020 EMD Millipore Corporation, an affiliate of Merck KGaA, Darmstadt, Germany. All rights reserved. Each product is sold with a limited warranty, which is provided with each purchase. Each product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes, nor is it intended for human use. EMD Millipore products may not be resold, modified for resale, or used to manufacture commercial products without written approval of EMD Millipore.

EMD Millipore Corporation | 28820 Single Oak Dr., Temecula, CA 92590

Technical Support NA +1-800-221-1975 | email: www.millipore.com/techservices | www.calbiochem.com

Technical Support All Other Countries - Contact Your Local Office

FOR RESEARCH USE ONLY.

Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purpose other than research is strictly prohibited.

Calbiochem and all other trademarks, unless specifically identified as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany