



Product Specification Sheet

VWR PIPETTE TIPS, RELOAD SYSTEM, RACKED

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Type	Capacity (µl)	Colour code	Packed	Pk
613-1075	Micro, graduated	0.1 - 10	Transparent	10 racks, 96 each	960
613-1076	Micro, extra long, graduated	0.1 - 10	Transparent	10 racks, 96 each	960
613-1074	Micro, Eppendorf type	0.5 - 20	Transparent	10 racks, 96 each	960
613-1072	Graduated, bevelled	1 - 200	Transparent	10 racks, 96 each	960
613-1073	Micro	1 - 200	Transparent	10 racks, 96 each	960
613-1071	Graduated, bevelled	5 - 300	Transparent	10 racks, 96 each	960
613-1070	Graduated	100 - 125	Transparent	10 racks, 96 each	960

Other Characteristics:

- Non sterile tips are free from DNA, DNase and RNase.
- Robust, traceable racks are easy to grab with their ergonomic clips
- Reload boxes including racks or empty reload boxes are also available and feature flip lock latch for one handed opening and closing, reversible lid with lip, flush fitting and click clip to stabilise the racks.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **October 2013**

