

Certificate of Conformity

VWR Catalogue Number	VWRA734-3330
Description	VWR®, Classic, PCR 96-well plate sealing membrane
Lot Number	20707/2526_01
Date of Manufacture (yyyy-mm-dd)	2025-01
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2026-05-07

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO ISO 13485:2016**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene (PP). Resins are USP Class-VI certified, FDA regulation CFR 21, and are FDA approved for food contact.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination. Test limit is 1.0×10^{-6} Kunitz Units for DNase & 1.0×10^{-9} Kunitz Units for RNase.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE) as defined in EC Directive 1069/2009/EC and EC Directive 142/2011/EC.
Cytotoxicity	This product is free of cytotoxicity effects.
Non-Pyrogenic Statement	Endotoxin levels are tested by Limulus Amoebocyte Lysate (LAL) Assay. Product is free from Endotoxin. The acceptance level for product is ≤ 0.01 EU/ ml (LAL Gel Clot Method).
Latex Statement	Product is free from latex.
BPA Statement	Bisphenols has not been intentionally incorporated in the manufacturing of the products.
DEHP Statement	DEHP has not been intentionally incorporated in the manufacturing of the products.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr^{6+}), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Bis(2-Ethylhexyl) Phthalate (DEHP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Ambient Room temperature.
Shelf Life	3 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.