

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® CLEANROOM COVERALL WITH HOOD STERILE
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2023-05-02

Catalogue Number	Colour	Clothing Size
114-9324	White	S
114-9325	White	M
114-9326	White	L
114-9327	White	XL
114-9328	White	XXL
114-9329	White	3XL
114-9330	White	4XL

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

	Test Method	MD	CD
Break Tensile properties(N)	ASTM D5034	> 60	> 35
Elongation rate [%]	ASTM D5034	> 30	> 30
Tear prosperities	ISO9073-4	> 30	> 30

Product Specifications

Material	60g SMS (Spunbonded Polypropylene fabric /Meltblown /Spunbonded Polypropylene fabric)
Sterilization	Products are Gamma irradiated as per ISO 11137:2006 and meets SAL 10 ⁻⁶ with a specified dose of 27.6 – 41.9 kGY.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	0- 50 °C, 45% - 65% humidity, Shady
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.