

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Magnetic stirring discs
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2024-10-21

Catalogue Number	Length [mm]	Diameter [mm]	Pk
442-0490	27	27	5
442-0491	37	37	5
442-0492	37	37	5

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Virgin Polytetrafluoroethylene (PTFE) Grade, strong Alnico V magnetic core. Raw material is in compliance with Directive FDA 21 CFR 177.1550 and complies with USP Class VI in accordance with USP <88> and <87>.
Temperature Resistance	–200 to +280 °C.
Sterilization	Non-Sterile.

ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP are not used in the manufacture of the raw material and are not expected to be present.
RoHS	Products are in compliance with EU Directive 2011/65/EU as amended by EU 2015/863.
REACH Statement	Product meet the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006.
Storage Conditions	Ambient temperature.
Shelf Life	Indefinite when stored under ideal conditions.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	10-2024