



VWR® 25mL Tissue Culture Flasks

Description: 25mL, Non-Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is not desired


Materials

Flask: GPPS (General polystyrene) Color: Clear

Cap: HDPE(High Density Polyethylene) Color: White

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Dimensions (Unit: mm)



Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 8 mL
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system




VWR® 25mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2783	10861-566	25mL Tissue Culture Flask, Non-Treated, w/plug seal cap, Sterilized	Yes	10	50	200
734-2784	10861-568	25mL Tissue Culture Flask, Non-Treated, w/vent cap, Sterilized	Yes	10	50	200



Dimensions (Unit: mm)



VWR® 50mL Tissue Culture Flasks

Description: 50mL, Non-Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is not desired

Materials

Flask: GPPS (General polystyrene) Color: Clear

Cap: HDPE(High Density Polyethylene) Color: White

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10® (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 17.5 mL
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE

R

VWR® 50mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2785	10861-570	50mL Tissue Culture Flask, Non-Treated, w/plug seal cap, Sterilized	Yes	10	50	200
734-2786	10861-572	50mL Tissue Culture Flask, Non-Treated, w/vent cap, Sterilized	Yes	10	50	200



VWR® 250mL Tissue Culture Flasks

Description: 250mL, Non-Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is not desired

Materials

Flask: GPPS (General polystyrene) Color: Clear


Cap: HDPE(High Density Polyethylene) Color: White

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 60 mL
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



STERILE R

VWR® 250mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2787	10861-574	250mL Tissue Culture Flask, Non-Treated, w/plug seal cap, Sterilized	Yes	5	25	100
734-2788	10861-576	250mL Tissue Culture Flask, Non-Treated, w/vent cap, Sterilized	Yes	5	25	100



VWR® 600mL Tissue Culture Flasks

Description: 600mL, Non-Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is not desired

Materials

Flask: GPPS (General polystyrene) Color: Clear


Cap: HDPE(High Density Polyethylene) Color: White

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 125 mL
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)




VWR® 600mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2789	10861-578	600mL Tissue Culture Flask, Non-Treated, w/plug seal cap, Sterilized	Yes	5	20	40
734-2790	10861-580	600mL Tissue Culture Flask, Non-Treated, w/vent cap, Sterilized	Yes	5	20	40



Dimensions (Unit: mm)



VWR® 600mL Extended Tissue Culture Flasks

Description: 600mL, extended, Non-Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is not desired

Materials

Flask: GPPS (General polystyrene) Color: Clear

Cap: HDPE(High Density Polyethylene) Color: White

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Extended height reduces the risk of overflowing media
- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 200 mL
- Meet the requirements of BSE/TSE
- Latex Statement: this product is latex free
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 600mL Extended Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/PK	Qty/Box	Qty/Case
734-3252	76211-292	600mL Tissue Culture Flask, extended, Non-Treated, w/plug seal cap, Sterilized	Yes	5	20	40
734-3253	76211-294	600mL Tissue Culture Flask, extended, Non-Treated, w/vent cap, Sterilized	Yes	5	20	40
	76308-724	Plug Seal Cap for 600mL Extended Tissue Culture Flask, Sterilized	Yes	5	50	200



VWR® 850mL Tissue Culture Flasks

Description: 850mL, Non-Treated, w/plug seal cap or w/vent cap, Sterilized

Purpose: For use in applications where cell attachment is not desired

Materials

Flask: GPPS (General polystyrene) Color: Clear


Cap: HDPE(High Density Polyethylene) Color: White

Membrane: 0.22µm Hydrophobic PP (Only use for Vent Cap)

Features

- Temperature: -20°C~+50°C
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- DNase/RNase free and Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Maximum working volume: 200 mL
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



STERILE R

VWR® 850mL Tissue Culture Flasks

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Bag	Qty/Pk	Qty/Case
734-2791	10861-582	850mL Tissue Culture Flask, Non-Treated, w/plug seal cap, Sterilized	Yes	3	9	18
734-2792	10861-584	850mL Tissue Culture Flask, Non-Treated, w/vent cap, Sterilized	Yes	3	9	18