




Sol - Vet® Butterfly Needles

- Removeable luer adapter allowing device to be used for blood draw or infusion
- Thin-walled needle minimizes patient discomfort
- Transparent design allows for clear visualization of blood flow
- Available with removeable multi-sample luer adapter for use with infusion or blood collection

***samples available upon request**

VWR Cat No.	Product Description	Box Qty	Case Qty
76582-506	Sol-Vet Butterfly Needle with Adapter 19G x 3/4 x 12"	50	1,000
76582-500	Sol-Vet Butterfly Needle with Adapter 21G x 3/4 x 12"	50	1,000
76582-502	Sol-Vet Butterfly Needle with Adapter 23G x 3/4 x 12"	50	1,000
76582-504	Sol-Vet Butterfly Needle with Adapter 25G x 3/4 x 12"	50	1,000