



MEDICATION DELIVERY
NORTH COVE FACILITY
Hwy 221 N PO Box 1390
Marion N.C. 28752
Telephone: (828) 756-4151
Fax: (828) 756-4821

Certificate of Analysis

Product: 0.9% Sodium Chloride Irrigation, USP

Lot #: C920272

Code: 2B7126

Manufacture/

Sterilization Date: 10/23/2013

Expiry: 10/2015

Chemical Testing per Specification: 03-15-19-061

Issue Date: 10/3/2008

TEST	LIMIT	RESULT
NaCl (g/L)	8.55 - 9.45 g/L	8.87 g/L
Sodium ID	Positive	Positive
Sodium ID- Flame	Positive	Positive
pH at 25 deg. C	4.5 - 7.0	5.6
Particle Analysis	NMT 25 $\geq$ 10 um/ml	NMT 25
Particle Analysis	NMT 3 $\geq$ 25 um/ml	NMT 3
Sterility	Pass Parametric Release	Pass

Limulus Testing per 13-01-V

RESULT	LIMIT	DISPOSITION	DATE PASSED
< 0.25 EU/ml	< 0.25 EU/ml	Pass	11/14/2013

This is to certify that this product was manufactured according to current GMP and fulfills the requirements of the Master Production Document.

Date Batch Released	Parametric Release Date	Quantity Released
11/14/2013	11/14/2013	NA

	Print Name	Signature/Date
Prepared by:	Sherrie Hopson	Sherrie Hopson 12/4/13
Approved by: Quality Manager or Designee	Keith Hawk	Keith Hawk 12/4/13