

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Flat Cleanroom Mops
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2025-06-25

Mop Head

Catalogue Number	Dimension [mm]	Sterile	Pk
129-0067	230	Non-Sterile	24
129-0068	356	Non-Sterile	24
129-0069	457	Non-Sterile	24
129-0070	230	Sterile	24
129-0071	356	Sterile	24

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	PolyTex, yellow foam.
Sterilization	Products are gamma irradiated and meets a minimum requirement of 10^{-6} SAL with a dose range of 25 – 40 kGy.(which id not applicable for non-sterile products)
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	This product is free of Cytotoxicity effects.
Non-Pyrogenic Statement	Products are pyrogen free.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP are not used in the manufacture of the raw material and are not expected to be present.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Product meet the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006. To the best of our knowledge product does not contain SVHC's above the regulatory threshold of 0.1%.
Storage Conditions	Ambient temperature recommended 15 °C – 30 °C. Do not store in direct sunlight. Store the product in a space that is protected from extreme temperature and humidity fluctuations and away from ultraviolet sources.
Shelf Life	Shelf life is not applicable non-sterile products and 5 years for sterile products.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	06-2025
2	1.1		
3	1.2		