

|                      |                                          |
|----------------------|------------------------------------------|
| Material             | 94006.1016                               |
| Material description | Ipasept                                  |
| Grade                | GAMMA-IRRADIATED 25 kGy 0,2 µm filtrated |
| CAS Number           | 67-63-0                                  |
| Molecular formula    | CH <sub>3</sub> CHOHCH <sub>3</sub>      |

| Characteristics                | Specifications          |
|--------------------------------|-------------------------|
| Assay (2-Propanol) (V/V)       | 68.0 - 72.0 %           |
| Assay (2-Propanol) (V/V) [USP] | Passes test USP         |
| Appearance                     | Clear colourless liquid |
| Acidity                        | Passes test USP         |
| Raw materials                  | Passes test Ph.Eur.     |
| Density (20/4)                 | 0.872 - 0.883           |
| Non-volatile matter            | ≤ 100 mg/l              |

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode  
VWR International bv; Geldenaaksebaan 464; BE-3001  
Leuven; Belgium

The irradiation process has been validated according to ISO 11137-2:2013, VDMAX 25 method.  
Starting from batch 14F a sterility assurance level of 10<sup>-6</sup> can be guaranteed on the entire content of the inner bag.

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.