




REVISION NO	REVISION DATE	DESCRIPTION	REVISED BY	APPROVED BY

SQUARE WIDE MOUTH BOTTLE SPECIFICATION

VWR-US CAT NO	VWR-EUROPE CAT NO	MATERIAL	COLOR	CAPACITY (ml)	APROX BRIM CAPACITY (ml)	CAP SIZE	HEIGHT WITH OUT CAP (mm)	BODY WIDTH (mm)	BODY LENGTH (mm)
414004-190	215-1904	HDPE	TRANSLUCENT	4000ml	4560ml	100-415	290	147	147
414004-189	-	PP	TRANSLUCENT	4000ml	4560ml	100-415	290	147	147

NOTE :

- 1) PP MATERIAL BOTTLES ARE AUTOCLABLE.
- 2) BOTTLE HAS INTEGRATED HANDLE DESIGN.

Closure Application Torque		
Closure Size	Minimum Torque Kg cm	Maximum Torque Kg cm
13-415	6	8
20-415	11	16
24-415	13	19
28-415	18	25
38-415	25	35
38-430	31	38
43-415	32	44
48-415	34	48
53-415	38	52
53T	43	60
63-415	46	64
83T	69	96
100-415	54	74

TITLE : WIDE MOUTH SQUARE BOTTLE

Material : AS PER TABLE

Rev : 0

SCALE : NOT TO SCALE



DRAWN BY	
DRAWN DATE	2/2/2009
DRAWING NO	

