

Certificate of Conformity

VWR Catalogue Number	See List
Description	Safety Wash Bottles with Integrated Safety Vent Cap
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2021-11-05

Safety wash Bottle, Wide Neck 250 mL

Catalogue Number	Imprint	Type
215-4294	Acetone	Natural bottle, wide neck, red cap
215-4295	Methylene dichloride	Natural bottle, wide neck, yellow cap
215-4296	Distilled water	Natural bottle, wide neck, natural cap
215-4297	Ethanol	Natural bottle, wide neck, white cap
215-4298	Isopropanol	Natural bottle, wide neck, yellow cap
215-4299	Methanol	Natural bottle, wide neck, green cap
215-4300	Methylethylketone	Natural bottle, wide neck, green cap
215-4301	Sodium hypochlorite	White bottle, wide neck, blue cap
215-4302	Toluene	Red bottle, wide neck, red cap
215-4303	Unprinted, with inscription label	Natural bottle, wide neck, white cap

Safety wash Bottle, Wide Neck 500 mL

Catalogue Number	Imprint	Type
215-4304	Acetone	Natural bottle, wide neck, red cap
215-4305	Methylene dichloride	Natural bottle, wide neck, yellow cap
215-4306	Distilled water	Natural bottle, wide neck, natural cap
215-4307	Ethanol	Natural bottle, wide neck, white cap
215-4308	Isopropanol	Natural bottle, wide neck, yellow cap
215-4309	Methanol	Natural bottle, wide neck, green cap

Catalogue Number	Imprint	Type
215-4310	Methylethylketone	Natural bottle, wide neck, green cap
215-4311	Sodium hypochlorite	White bottle, wide neck, blue cap
215-4312	Toluene	Red bottle, wide neck, red cap
215-4313	Unprinted, with inscription label	Natural bottle, wide neck, white cap

Safety wash Bottle, Wide Neck 1000 mL

Catalogue Number	Imprint	Type
215-4314	Acetone	Natural bottle, wide neck, red cap
215-4315	Methylene dichloride	Natural bottle, wide neck, yellow cap
215-4316	Distilled water	Natural bottle, wide neck, natural cap
215-4317	Ethanol	Natural bottle, wide neck, white cap
215-4318	Isopropanol	Natural bottle, wide neck, yellow cap
215-4319	Methanol	Natural bottle, wide neck, green cap
215-4320	Methylethylketone	Natural bottle, wide neck, green cap
215-4321	Sodium hypochlorite	White bottle, wide neck, blue cap
215-4322	Toluene	Red bottle, wide neck, red cap
215-4323	Unprinted, with inscription label	Natural bottle, wide neck, white cap

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** and **ISO 13485:2016** standards. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-05	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	NA
Dimensional	QSF-QCIR-02	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

Material	Bottle: Low Density Polyethylene (LDPE), 100 % Virgin, FDA 21CFR Compliance. Closure and Draw tube: Polypropylene (PP), 100 % Virgin, FDA 21CFR Compliance.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.

Storage Conditions Store at Room temperature.

Shelf Life Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.