

Certificate of Conformity

VWR Catalogue Number	See list
Description	Slider Zip bags
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-11-21

Slider zip bags, transparent, with blue slider, without printed area

Catalogue Number	Thickness [μm]	Dimension [mm]
129-1198	50	250 × 350
129-1199	50	300 × 400

Slider zip bags, transparent, with blue slider, without printed area, with bottom gusset

Catalogue Number	Thickness [μm]	Dimension [mm]
129-1202	50	250 × 350
129-1203	50	300 × 400

Slider zip bags, transparent, with blue slider, with write-on area

Catalogue Number	Thickness [μm]	Dimension [mm]
129-1200	50	250 × 350
129-1201	50	300 × 400

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	LDPE 260GG + FD21SH
Heat resistant	-10 to +65 °C.
Tensile Strength	> 21 MP/10mm.
Weld Strength	> 10N/15mm.
Elongation	> 550%.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.

Storage Conditions Room Temperature.

Shelf Life No shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	11-2022