

Certificate of Conformity

VWR Catalogue Number	See list
Description	Disposal bags, Autoclavable
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-03-09

Without 'Biohazard' printing

Catalogue Number	Description	Capacity [L]	Length x Width [mm]
129-0578	Autoclavable bags, clear	2	300 × 200
129-0579	Autoclavable bags, clear	3	400 × 250
129-0580	Autoclavable bags, clear	12	660 × 305
129-0582	Autoclavable bags, clear	14	610 × 420
129-0581	Autoclavable bags, clear	16	780 × 400
129-0583	Autoclavable bags, clear	25	770 × 550
129-0584	Autoclavable bags, clear	32	780 × 600
129-0585	Autoclavable bags, clear	40	1100 × 700
129-0586	Autoclavable bags, clear	60	1100 × 900

With 'Biohazard' printing

Catalogue Number	Description	Capacity [L]	Length x Width [mm]
129-0587	Autoclavable bags, clear, with 'Biohazard' printing	12	660 × 305
129-0588	Autoclavable bags, clear, with 'Biohazard' printing	14	610 × 420
129-0589	Autoclavable bags, clear, with 'Biohazard' printing	32	780 × 600

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene HP528N (Virgin).
Heat resistant	Up to 145 °C.
Tensile Strength	> 25 MP/10mm.
Weld Strength	> 12N/15mm.
Elongation	> 600%.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.

RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room Temperature.
Shelf Life	No shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.