

PRODUCT DATA SHEET

ORIG.PERFUSOR SYR. 50ML LL (8728844F)



Attribute Name	Attribute Value
Brand name	Perfusor®
Product name (HC)	Original Perfusor® Syringe
Product type	Syringe
Application	Infusion pump
Construction of syringe	3-piece
Style of syringe cone	Luer Lock
Position of syringe cone	centric
Nominal capacity	50.0 ml
Total graduated capacity	50.0 ml
Max. usable capacity	50.00 ml
Scale value	1.0000 ml
Light-protected	No

Attribute Name	Attribute Value
Color of light protection	-
Catalogue attribute	without aspiration needle
Type of packaging	Blister
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Product contains BPA	No
Lubricants	Silicone oil
Reference norm	ISO 11135-1
Sterility	sterile
Type of Sterilization	EO
Disposable / Reusable	disposable