

Material No.: T007-09
Revision No.: 0

Product Specification

Test	Specification
Assay ($\text{CH}_3(\text{CH}_2)_3\text{CH}_3$) (by GC, corrected for water)	$\geq 98 \%$
Color (APHA)	≤ 10
Residue after Evaporation	$\leq 0.0020 \%$
Sulfur Compounds (as S)	$\leq 0.005 \%$
Water (by KF, coulometric)	$\leq 0.02 \%$
Boiling Range (95%)Lower	34 – 36 °C
Boiling Range (95%)Upper	34 – 40 °C

For Laboratory, Research, or Manufacturing Use
Meets Reagent Specifications for testing USP/NF monographs

Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1