

Product Data Sheet

Commercial grade glass hydrometer


Reference 3020FG140/20-qp

For general purpose

Density hydrometer (g/ml)

Reading below meniscus

Not available for Calibration
Certificate



Specifications

- Range 1400-1600 g/ml
- Division: 0.002 g/ml
- Accuracy: +/- 0.004 g/ml
- Calibration temperature : 20°C
- Surface tension : 75 mN/m
- ± 300 mm
- In individual box

Alla France

Z.I. du Bompas - 13, rue du Bompas - Chemillé - 49120 Chemillé-en-Anjou - FRANCE

Tél. + 33 (0)2 41 30 55 08 - Fax +33 (0)2 41 30 34 67

www.allafrance.com

