



VWR® Aspirating Pipettes

Description: Aspirating Pipettes

Purpose: Transfer liquid


Materials: 100% high clarity polystyrene

Color: Clear

Features

- Uses: under -20 ~ 50 °c
- Sterilization: sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- BSE/TSE Statement: Meet the requirements of BSE/TSE.
- Products are manufactured under the ISO 13485 standard.
- Shelf Life: 4 years after month of production

VWR® Aspirating pipettes



EU Part No.	Volume	Total Length (mm) ± 2mm	Tip		Mouth		Sterile SAL 10 ⁻⁶	individually Pack	Packing
	(ml)		Outside Dia.	interior Dia.	Outside Dia.	interior Dia.			
612-5883	1	278	2.30	1.50	4.80	3.00	YES	YES	50pcs/bag, 200pcs/box
612-5884	2	278	2.50	2.00	6.28	4.26	YES	YES	50pcs/bag, 200pcs/box
612-5885	5	347	3.30	1.53	8.00	6.00	YES	YES	50pcs/bag, 200pcs/box
612-5886	10	343	3.40	1.50	8.00	5.43	YES	YES	50pcs/bag, 200pcs/box
612-5887	25	312	4.75	2.35	8.00	5.55	YES	YES	50pcs/bag, 200pcs/box

