

BAKERBOND spe, SDB

BAKERBOND spe

Disposable Extraction Columns; Jumbo pack



Material No: 7519-22
Effective from: 2023/11/08

Product Specification

Test	Specification
Styrene Divinylbenzene (SDB-1):	
Breakthrough Antrazine*	$\geq 1 \text{ L}$
Breakthrough Desisopropylatrazine*	$\geq 1 \text{ L}$
Specific Surface Area	$\geq 950 \text{ m}^2/\text{g}$
Porosity, Median Pore Diameter	50 – 70 Å
Particle Size – 50% Volume	65 – 95 μm
Capacity (mg caffeine / g sorbent)	≥ 250
Recovery Atrazine*	$\geq 97 \%$
Recovery Desisopropylatrazine*	$\geq 75 \%$

For Laboratory, Research or Manufacturing use.
3mL Solid Phase Extraction Columns
200 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000