



CERTIFICATE OF ANALYSIS

BAC 20WN[®]
Gel

Lot # 13420A

Date of Manufacture May 2020

Recommended Retest date is two years from the date of manufacture

SIZE	4oz (120mL)	16oz (480mL)	1gal (3.8L)
DECON #	7104	7106	7101

TEST	SPECIFICATION	UNIT	TEST RESULT
Appearance	Clear to slightly hazy viscous liquid	N/A	Pass
pH	7 - 8	pH	7.6
Viscosity @ 21- 23° C	8,000 to 12,000	cps	11,447

Reviewed and Approved by:

Date: May 2020

A handwritten signature in blue ink, appearing to read 'Karen Wend', is written over a light blue horizontal line.

Karen Wend
Manager, Quality and Technical Services

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.