



Specification

1.07224.1000 Peptone from meat peptic digested, granulated, for microbiology

Specification		
Appearance	Lightly brownishyellow granulate.	
Identity	passes test	
Total nitrogen (N) (Kjeldahl)	10.0 - 16.0	%
Amino nitrogen (as N)	≥ 1.5	%
pH-value (5 %; water)	6.5 - 7.5	
Nitrite (NO ₂)	passes test	
Sulfated ash (800 °C)	≤ 15.0	%
Loss on Drying (105 °C)	≤ 6.0	%
Suitability for microbiology	passes test	

Dr. Lukas Mechler

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.