

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Jars, Wide Mouth, with Screw Cap
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2026-06-01

Catalogue Number	Capacity [mL]	Height [mm]	Ø cap [mm]	Pack
216-0164	15	34	34	1
216-0165	30	62	34	1
216-0166	60	46	60	1
216-0167	120	62	72	1
216-0168	240	100	72	1
216-0169	360	95	90	1
216-0170	480	125	90	1
216-0171	1000	160	110	1
216-0172	2200	260	134	1

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Virgin Polytetrafluoroethylene (PTFE) Grade PTFE 807N X.
Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI plastic (121 °C).

Sterilization Non-Sterile.

ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	Product does not contain any of the restricted substance above the maximum concentration values as listed in Annex II of RoHS Directive 2011/65/EU as amended by EU 2015/863.
REACH Statement	Product does not contain any Substance of very high concern (SVHC) at a concentration of $\geq 0.1\%$ (w/w) that are on European Union (EU) REACH candidate list as published by the European Chemical Agency (ECHA).
Storage Conditions	Ambient temperature.
Shelf Life	Shelf life is generally considered to be indefinite, within the confines of normal use and storage conditions.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	01-2025
2	2.0	Modified Quality system compliance and disclaimer	08-2025
3	2.1	Modified USP Class VI information	06-2026