



Certificate of Conformity

VWR® BASIC MOBCAP DOUBLE ELASTIC DETECTABLE



CHARACTERISTICS:

- Made of non-woven SPP fabric
- Multiple colours option from white, blue, green to red.
- Nature Rubber Latex Free
- Double elastic
- Detectable

PHYSICAL PROPERTIES

Item		Test method	Standard
Weight(g/m ²)		ISO3801	10 ± 1
Break Tensile Properties (N/5cm)	MD	ISO13934-1	>15
	CD	ISO13934-1	>8
Elongation Rate(%)	MD	ISO13934-1	>30
	CD	ISO13934-1	>30

PRODUCTS AND DIMENSION SPECIFICATION



Unit:CM

ECN		113-8895	113-8896	113-8897	113-8898	Tolerance
Colour		White	Blue	Green	Red	
Dimension	A	21.5	21.5	21.5	21.5	± 2
	A'	53	53	53	53	± 1
	B	26	26	26	26	± 1

Remark: A: Relax, A': Stretch



Certificate of Conformity

PACKING:

Package/Quantity	100pcs/bag, 10bags/case
Best Storage Conditions	Temp : 0- 50°C, Humidity 45-85% without sunshine exposure
Country of Origin	China

QUALITY ASSURANCE CRITERIA:

1. Dimension and Function - Critical (AQL1.5, S-2)

- Product dimension out of spec.
- Fabric characteristics out of spec
- Sonic welding strength less than 15N/5cm

2. Visual Inspection/ Attributes - Major (AQL 2.5, L-II)

- Foreign Material
- Dirt
- Thread ends loose or too long(over 1cm)
- Holes/Tears
- Foul odor
- Discolor
- Wrong collar tag
- Wrong stitch
- Cloth surface applied wrongly
- Less than 3 stitches per centimeter

3. Visual Inspection/Attributes - Minor (AQL 4.0, L-II)

- Sewing sundries
- Elastic band end exposed
- Sewing end falling-off
- Jumping Stitch
- Missing stitch
- Unwanted folds
- Elastic tightness non-consistent

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**