

Material	271402T
Material description	Ammonia 25%
Grade	GPR RECTAPUR

CAS Number	1336-21-6
Molecular formula	H3N
Molecular mass	17.03

Characteristics	Specifications
Assay	≥ 25 %
Non-volatile residue	≤ 100 ppm
Cl (Chloride)	≤ 50 ppm
SO ₄ (Sulphate)	≤ 50 ppm
Fe (Iron)	≤ 5 ppm
Pb (Lead)	≤ 5 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.