



Certificate of Conformity


VWR Plain Microscope Slide

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Description	Thickness	LxW	Pk
631-1160	Ground edges, clipped corner (45°) slides	1.0 (-0.05/+0.1) mm	76×25 mm	2.500



Purpose: Used for hematology, microbiology

Material: Standard Glass

Features: The smooth edges provide the safety of the operation

Temperature: Normal laboratory temperature and humidity

Uses: For the human body blood or fluid test

Shelf life: 60 Months

This document has been produced electronically and is valid without a signature

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2019**