

Material	25656.298
Material description	4-Methyl-2-pentanone
Grade	GPR RECTAPUR
CAS Number	108-10-1
Molecular formula	$(CH_3)_2CHCH_2COCH_3$
Molecular mass	100.16

Characteristics	Specifications
Assay	$\geq 98 \%$
IR Spectrum	Passes test
Boiling point	115 - 117 °C
Density (20/4)	0.790 - 0.810
Evaporation residue	$\leq 0.02 \%$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.