

Certificate of Conformity

VWR Disposable Cuvette, Macro

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 634-0675

Lot No: 0000002414

Type	Range	Window material	Light path	Volume	Pk
Macro	340 - 900 nm	PS	10 mm	2.5 – 4.5 ml	100

Each case of disposable cuvettes contains cuvettes of one cavity of the injection mold to ensure minimal variations of extinction values. The number of cavities is displayed on each cuvette.

Standard deviation at 360nm $\leq \pm 0.005$ extinction units.

This document has been produced electronically and is valid without a signature

Date of issue: **June 2021**